

# Broncoflex™

Single-use Bronchoscope



## Easy BAL Clip™

**TSC** The Surgical  
Company

— ENDOVISION —

Specialization through Innovation

# Easy BAL Clip™

## EASY SAMPLING SOLUTION

Easy BAL Clip™ is a single-use sampling solution for broncho-alveolar lavage procedure. Attached to the Broncoflex, it can be connected to multiple sizes of specimen containers.



## COMPATIBILITY WITH SPECIMEN CONTAINER

Easy BAL Clip™ allows connection to multiple sizes of specimen containers. No limit in the sample size or number of sampling.



## SINGLE OPERATOR SAMPLING SOLUTION

A physician can solo perform sampling, and switch to another container as needed.



## IMPROVED SUCTION LINE MANAGEMENT

Top handle placement keeps specimen container upright to prevent sample loss.

## TECHNICAL SPECIFICATIONS

MATERIAL	ABS-PC
LENGTH	90mm
WIDTH	75mm
HEIGHT	30mm
WEIGHT	15.5g
PART NUMBER	41010000

The Surgical Company Endovision  
Axess Vision Technology  
Zone de la Liodière - 6 rue de la Flottière  
37300 Joué-lès-Tours (France)  
Tel.: + 33 (0)2 47 34 32 90

<https://tsc-group.com/endovision>  
[www.broncoflex.com](http://www.broncoflex.com)

