



The Outdoor Cabinet for **All AEDs**

The AIVIA V3 PROTECTION is an eco-designed outdoor wall-mounted cabinet, compatible with all brands and models of AEDs (Automated External Defibrillators).

Designed to meet the new EU-MDR regulations, it ensures protection, compliance, and proper operation of AEDs 24/7, in all weather conditions.

Medical Device / EU-MDR – Temperature Control



Medical Device EU-MDR *

An AED cabinet is considered a medical device if it directly affects the AED's ability to function properly.

* European Medical Device Regulation

Technical Specifications



DURABILITY & RELIABILITY

- IK10 high impact resistance
- Sturdy materials: ABS-PC frame and UL94-V0 flame-retardant, UV-protected polycarbonate door



CERTIFICATIONS

- Fully compliant with the EU Medical Device Regulation (EU-MDR 2017/745)
- Certified under IEC 60601-1, the international standard for medical electrical equipment safety and performance



DIAGNOSTICS & MAINTENANCE

- Weekly self-tests, recorded and accessible for maintenance
- Event traceability via NFC-enabled smartphone access



VISUAL & AUDIBLE ALERTS

- Bright LED indicators and adjustable 90 dB buzzer
- Adaptive lighting: constant or light-sensor triggered



CLIMATE CONTROL – AUTOMATED

- Efficient heating system for outdoor temperatures as low as -20°C
- Built-in ventilation system
- IP54 weatherproof rating



ECO-DESIGN

- Made from recyclable green ABS and polycarbonate materials
- Optimized palletization to reduce environmental impact
- PYRESCOM is ISO 14001 certified and EcoVadis rated for its CSR commitments



PYRES City

PYRES Health

PYRES Entreprise



About AIVIA Cabinets

Since 2006, AIVIA has been protecting AEDs worldwide. Our cabinets suit all AED models and can be enhanced with management and monitoring tools to ease fleet maintenance and ensure device readiness.



700 000+
AIVIA cabinets
installed



Present in **50+**
countries



Over **19 years**
of experience



Système de
management
ISO 9001:2015
ISO 14001:2015
www.tuv.com
ID 9105060867



PYRES.COM

Contact us for more information !

+33 4 68 68 39 68

aivia@pyres.com

<https://pyres.com/en/solutions/aivia/>