



Technical Data Sheet

AIVIA V3 200



ELECTRICAL CONSUMPTION



- Minimum consumption: 100 mA
- Maximum consumption: 2700 mA (with heating active)

TEMPERATURE



Operating Temperature :

Minimum (static): -20°C

Maximum: +55°C

Temperature alerts and controls :

Low temperature warning : 5°C ± 2°C

High temperature warning : 45°C + -2°C

Heating activation : 10°C ± 2°C

Fan activation : 30°C ± 2°C

EUMDR COMPLIANCE



The AIVIA **complies with EU Regulation 2017/745 (EU-MDR)** on medical devices and meets the essential requirements of Annex I of this regulation.

The AIVIA is a **Class I device**.

- IEC 60601-1: 2005,AMD1: 2012 AMD2:2020
- IEC 62366-1: 2020
- IEC 60601-1-2: 2014 2015 / AMD1 2020
- IEC 62311: 2019
- IEC 62304: 2006, AMD1: 2015
- IEC 60601-1-11: 2015, A1: 2020
- EN 300 330 V2.1.1

HUMIDITY



Maximum relative humidity : 90% without condensation

ATMOSPHERIC PRESSURE



700 à 1 060 hPa

SOUND POWER



90 db / 1 meter