



PROVIDING INSIGHTS TO
IMPROVE OUTCOMES

MOVING CARE **FORWARD**


RescueNet® CaseReview

Seamlessly manage your patient case data



IMMEDIATE POST-CALL DATA ACCESS

Quickly transmit essential clinical data, including CPR quality, ECG, and EtCO₂, directly from an X Series® Advanced monitor/defibrillator or ZOLL AED 3® BLS defibrillator. The push of a button ensures that key metrics are delivered to QA/QI personnel – no cards or cables are required.



TRAINING

Establish quality standards for training purposes. When used with adult or pediatric CPR electrodes, CaseReview makes it easy to set training goals to improve your team's performance. Target your training to address identified areas of improvement and compare training versus actual results to drive improved outcomes.





SAME-DAY DEBRIEFING

Review all the details and data of a sudden cardiac arrest (SCA) or respiratory distress event while the information is still fresh. Help your team learn what was successful and what needs improvement. The playback tab displays clinical information from your ZOLL® device, including all vital signs, clinical waveforms, and CPR quality, as well as ventilation quality.



AutoPulse® and X Series Advanced:


Two devices, one report

Merging data from multiple devices can be complicated and time-consuming. AutoPulse Resuscitation System and X Series Advanced integrate data smoothly without complex data transfers, data cards, and cables. One patient – one report.



SYSTEM-WIDE QI

Identify trends and areas for improvement by viewing aggregate data. Improving outcomes requires a big-picture view of performance data. CaseReview saves valuable time and resources by quickly aggregating critical data such as peri-shock pause times, CPR quality, manual ventilation quality, and more, displaying it in colorful, intuitive graphics.



RescueNet CaseReview: For a complete view of the patient journey

To provide the highest quality care every step of the way, your team needs complete information on the patient journey — from initial response to post-event data. CaseReview helps you achieve the best outcomes in SCA and respiratory distress calls by tracking and managing case data that can inform and improve patient outcomes and operational efficiencies. Whether collecting data for a hot debrief or identifying opportunities for system wide performance improvements during a cold debrief, CaseReview helps your team close that loop.

ZOLL MEDICAL CORPORATION

269 Mill Road | Chelmsford, MA 01824 | 978-421-9655 | 800-804-4356 | zoll.com

Copyright © 2021 ZOLL Medical Corporation. All rights reserved. RescueNet, X Series Advanced, AutoPulse, ZOLL AED 3, and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries. All other trademarks are the property of their respective owners.

Printed in U.S.A.
MCN EP 2104 0373

For subsidiary addresses and fax numbers, as well as other global locations, please go to zoll.com/contacts.

ZOLL®