




Intravascular Catheter Kit

for Thermogard™ Temperature Management Systems



- Streamlined for maximum efficiency and simplicity
- Reliable, precise catheter performance
- Dedicated kit with components essential for intravascular temperature control catheter use; accessories for venous access not included.
- **Easy to use with just 4 components:**
 - Catheter
 - Guidewire and dispenser
 - Introducer
 - Suture tab and clip
- Catheter options available for femoral or IJ/subclavian access, range of cooling powers



Contact your ZOLL® representative for more information or to learn which intravascular temperature management solution best supports your practice.

CATALOG #	DESCRIPTION
8700-0782-14	Icy® Temperature Control Catheter
8700-0783-14	Quattro® Temperature Control Catheter
8700-0793-14	Solex 7™ Temperature Control Catheter